


Product Name: Recombinant Human DBT Protein
Catalog #: PEH9118



Summary

Name	ODB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11182
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DIGEGIREVT VKEWYVKEGD TVSQFDSICE VQSDKASVTI TSTRYDGVIKK LYYNLDDIAY VGKPLVDIET EALKDSEEDV VETPAVSHDE HTHQEIKGRK TLATPAVRRRL AMENNIKLSE VVGSGKDGRI LKEDILNYLE KQTGAILPPS PKVEIMPPPP KPKDMTPVIL VSKPPVFTGK DKTEPIKGFQ KAMVKTMSAA LKIPHFGYCD EIDLTELVKL REELKPIAFA RGIKLSFMPF FLKAASLGLL QFPILNASVD

image



Background

Alternative Names	Lipoamide acyltransferase component of branched-chain alpha-keto acid dehydrogenase complex, mitochondrial (EC 2.3.1.168) (52 kDa mitochondrial autoantigen of primary biliary cirrhosis) (Branched chain 2-oxo-acid
--------------------------	--

Product Name: Recombinant Human DBT Protein
Catalog #: PEH9118



dehydrogenase complex component E2) (BCOADC-E2) (Branched-chain alpha-keto acid dehydrogenase complex component E2) (BCKAD-E2) (BCKADE2) (Dihydrolipoamide acetyltransferase component of branched-chain alpha-keto acid dehydrogenase complex) (Dihydrolipoamide branched chain transacylase) (Dihydrolipoyllysine-residue (2-methylpropanoyl)transferase)

Note

For Research Use Only , Not for Diagnostic Use