


Product Name: Recombinant Human RCC1 Protein
Catalog #: PEH9123



Summary

Name	RCC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P18754
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPKRIAKRRS PPADAIPKSK KVKVSHRSHS TEPGLVLTG QGDVGQLGLG ENVMERKKPA LVSIPEDVVQ AEAGGMHTVC LKSGQVYSF GCNDEGALGR DTSVEGSEMV PGKVELQEKV VQVSAGDSHT AALTDDGRVF LWGSFRDNNG VIGLLEPMKK SMVPVQVQLD VPVVKVASGN DHLVMLTADG DLYTLGCGEQ GQLGRVPELF ANRGGRQGLE RLLVPKCVML KSRGSRGHVR FQDAFCGAYF TFAISHEGHV YGFGLSNYHQ LGTPGTESCF IPQNLTSFKN STKSWVGFSG GQHHTVCMDS EGKAYSLGRA EYGRGLGEG AEEKSIPTLI SRLPAVSSVA CGASVGYAVT KDGRVFAWGM GTNYQLGTGQ DEDAWSPVEM MGKQLENRVV LSVSSGGQHT VLLVKDKEQS

image



Background

Product Name: Recombinant Human RCC1 Protein
Catalog #: PEH9123



Alternative Names

Regulator of chromosome condensation (Cell cycle regulatory protein)
(Chromosome condensation protein 1)

Note

For Research Use Only , Not for Diagnostic Use