


Product Name: Recombinant Human SMARCA1 Protein
Catalog #: PEH9125



Summary

Name	SMCA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P28370
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FFPPRLFELL EKEILYYRKT IGYKVPRNPD IPNPALAQRE EQKKIDGAEP LTPEETEEKE KLLTQGFTNW TKRDFNQFIK ANEKYGRDDI DNIAREVEGK SPEEVMEYSA VFWERCNELQ

image



Background

Alternative Names	Probable global transcription activator SNF2L1 (EC 3.6.4.-) (ATP-dependent helicase SMARCA1) (Nucleosome-remodeling factor subunit SNF2L) (SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 1)
--------------------------	---

Note



Product Name: Recombinant Human SMARCA1 Protein
Catalog #: PEH9125

For Research Use Only , Not for Diagnostic Use