


Product Name: Recombinant Human NUP62 Protein
Catalog #: PEH9128



Summary

Name	NUP62_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P37198
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAAASSAMTY AQLSLINKW SLELEDQERH FLQQATQVNA WDRTLIENTGE KITSLHREVE KVKLDQKRLD QELDFILSQQ KELEDLLSPL EELVKEQSGT IYLQHADEER EKTYKLAENI DAQLKRMAQD LKDIIIEHLNT SGAPADTSDP LQQICKILNA HMDSLQWIDQ NSALLQRKVE EVTKVCEGRR

image



Background

Alternative Names Nuclear pore glycoprotein p62 (62 kDa nucleoporin) (Nucleoporin Nup62)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NUP62 Protein
Catalog #: PEH9128

