

Product Name: Recombinant Human MLH1 Protein
Catalog #: PEH9130



Summary

Name	MLH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P40692
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKLDAFLQPL SKPLSSQPQA IVTEDKTDIS SGRARQQDEE MLELPAPAEV AAKNQSLEGD TTKGTSEMSE KRGPTSSNPR KRHREDS DVE MVEDDSRKEM TAACTPRRRI INLTSVLSLQ EEINEQGHEV LREMLHNHSF VGCVNPQWAL AQHQTKLYLL NTTKLSEELF YQILIYDFAN FGVRLRLSEPA PLFDLAMLAL DSPESGWTEE DGPKEGLAEY IVEFLKKKAE MLADYFSLEI DEEGNLIGLP LLIDNYVPPL EGLPIFILRL ATEVNWDEEK ECFESLSKEC

image



Background

Alternative Names DNA mismatch repair protein Mlh1 (MutL protein homolog 1)

Product Name: Recombinant Human MLH1 Protein
Catalog #: PEH9130



Note

For Research Use Only , Not for Diagnostic Use