


Product Name: Recombinant Human ATRX Protein
Catalog #: PEH9132



Summary

Name	ATRX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P46100
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKKKRRRIKV QEDSSSENKS NSEEEEEKE EEEEEEEEE EEEEDENDDS KSPGKGRKKI RKILKDDKLR TETQNALKEE EERRKRIAER EREREKLREV

image



Background

Alternative Names	Transcriptional regulator ATRX (EC 3.6.4.12) (ATP-dependent helicase ATRX) (X-linked helicase II) (X-linked nuclear protein) (XNP) (Znf-HX)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use