


Product Name: Recombinant Human ALDH1A3 Protein
Catalog #: PEH9133



Summary

Name	AL1A3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P47895
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GQPDRKPPAL PRPIRNLEVK FTKIFINNEW HESKSGKKFA TCNPSTREQI CEVEEGDKPD VDKAVEAAQV AFQRGSPWRR LDALSRGRLL HQLADLVERD RATLAALETM DTGKPFLHAF FIDLEGCI RT LRYFAGWADK IQGKTIPTDD NVVCFTRHEP IGVCGAITPW NFPLLMLVWK LAPALCCGNT MVLKPAEQTP LTALYLGLSI KEAGFPPGVV NIVPGFGPTV GAAISSHPQI NKIAFTGST VGKLVKEAAS RSNLKRVTLE LGGKNPCIVC ADADLDLAVE CAHQGVFFNQ GQCCTAASRV FVEEQVYSEF VRRSVEYAKK RPVGDPFDVK TEQGPQIDQK QFDKILELIE SGKKEGAKLE CGGSAMEDKG LFIKPTVFSE VTDNMRIAKE EIFGPVQPIL KFKSIEEVIK RANSTDYGLT AAVFTKNLDK ALKLASALE GTVWINCYNA LYAQAPFGGF KMSGNGRELG EYALAEYTEV KTVTIKLGDK

image



Background

Product Name: Recombinant Human ALDH1A3 Protein
Catalog #: PEH9133



Alternative Names

Aldehyde dehydrogenase family 1 member A3 (EC 1.2.1.36) (Aldehyde dehydrogenase 6) (Retinaldehyde dehydrogenase 3) (RALDH-3) (RaLDH3)

Note

For Research Use Only , Not for Diagnostic Use