


Product Name: Recombinant Human STAR Protein
Catalog #: PEH9135



Summary

Name	STAR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P49675
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EETLYSDQEL AYLQQGEEAM QKALGILSNQ EGWKKESQQD NGDKVMSKVV PDVGKVFRL E VVVDQPMERL YEELVERMEA MGEWNPVKE IKVLQKIGKD TFITHELAAE AAGNLVGPRD FVSVRCARR GSTCVLAGMA TDFGNMPEQK GVIRAEHGPT CMVLHPLAGS PSKTKLTWLL SIDLKGWLPK SIINQVLSQT QVDFANHLRK RLESHPASEA RC

image



Background

Alternative Names	Steroidogenic acute regulatory protein, mitochondrial (StAR) (START domain-containing protein 1) (StARD1)
--------------------------	---

Note

Product Name: Recombinant Human STAR Protein
Catalog #: PEH9135



For Research Use Only , Not for Diagnostic Use