


Product Name: Recombinant Human MSH6 Protein
Catalog #: PEH9138



Summary

Name	MSH6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P52701
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GTITYVTDKSE EDNEIESEEE VQPKTQGSRR SSRQIKRRV ISDES DIGG SDVEFKPDTK EEGSSDEISS GVG DSESEGL NSPVKVARKR KRMVTGNGSL KRKSSRKETP SATKQATSIS

image



Background

Alternative Names	DNA mismatch repair protein Msh6 (hMSH6) (G/T mismatch-binding protein) (GTBP) (GTMBP) (MutS protein homolog 6) (MutS-alpha 160 kDa subunit) (p160)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MSH6 Protein
Catalog #: PEH9138

