

Product Name: Recombinant Human OTUD4 Protein
Catalog #: PEH9147



Summary

Name	OTUD4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q01804
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSKIVMELDT LEVAEDNSE ISDEDDSC SKTAAAAADV NGFKPLSGNE QLKNNGNSTS LPLSRKVLKS LNPAVYRNVE YEIWLKSKQA QQKRDYSIAA GLQYEVGDKC

image



Background

Alternative Names	OTU domain-containing protein 4 (EC 3.4.19.12) (HIV-1-induced protein HIN-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use