


Product Name: Recombinant Human DEUP1 Protein
Catalog #: PEH9148



Summary

Name	DEUP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05D60
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HEKELNKIRS QLQQVEEYHN SEQERMNEI SDLTEELHQQ EITIATVTKK AALLEKQLKM ELEIKEKMLA KQKVSDMKYK AVRTENTHLK GMMGDLDPGE YMSMDFTNRE QSRHTSINKL QYENERLRND LAKLHVNGKS TWTNQNTYEE TGRYAYQSQI KVEQNEERLS HDCEPNRSTM PPLPPSTFQA KEMTSPLVSD DDVFPLSPPD MSFPASLAAQ HFLLEEEKRA KELEKLLNTH

image



Background

Alternative Names Deuterosome assembly protein 1 (Coiled-coil domain-containing protein 67)

Note

Product Name: Recombinant Human DEUP1 Protein
Catalog #: PEH9148



For Research Use Only , Not for Diagnostic Use