


Product Name: Recombinant Human EXOSC9 Protein
Catalog #: PEH9149



Summary

Name	EXOS9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q06265
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RNSKIDTES LCVVAGEK V W QIRVDLHLLN HDGNIIDAAS IAAIVALCHF RRPDVSVQGD EVTLTYPEER DPVPLSIHHM PICVSFAFFQ QGTYLLVDPN EREERVMDGL LVIAMNKHRE ICTIQSSGGI MLLKDQVLRC SKIAGVKVAE ITELILKALE NDQKVRKEGG KFGFAESIAN QRITAFKMEK APIDTSDVEE KAEEIIAAE PPSEVVSTPV LWTPGTAQIG EGVENSWGDL EDSEKEDDEG GGDQAII LDG IKMDTGVEVS DIGSQDAPII LSDSEEEEMI ILEPDKNPKK IRTQTTS AKQ EKAPSKKPVK

image



Background

Alternative Names	Exosome complex component RRP45 (Autoantigen PM/Scl 1) (Exosome component 9) (P75 polymyositis-scleroderma overlap syndrome-associated
--------------------------	--

Product Name: Recombinant Human EXOSC9 Protein
Catalog #: PEH9149



autoantigen) (Polymyositis/scleroderma autoantigen 1)
(Polymyositis/scleroderma autoantigen 75 kDa) (PM/Sci-75)

Note

For Research Use Only , Not for Diagnostic Use