


Product Name: Recombinant Human DNAH14 Protein
Catalog #: PEH9151



Summary

Name	DYH14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q0VDD8
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSLDEFCEEQ LQQATQALKQ LEDIRNKAIS EMKSTFLKVA EKNEIKEYFE SKLSEDDTTH FKLPHYRRL ETFFKFVMLV DYIFQELIRQ LMNTAVTLLL ELFNGSAGMP FSVEKKENEL IRTFKDNSFP TGKTTNDCEE LVDNSKLHAI SVQKSEVKTD TDINEILNSV EVGKDLRKTY APIFEVNLCL RIPAESDSSE NSKENFHESD QCPEECVMFE DEMSENKDNC VKKHSSEELL PKAKKSKEIS YNLEDIISAT ITPLCQDPQL SIFIDLVSIM DLPNKTGSII HYKEQTRWPD CHILFETDPA

image



Background

Alternative Names	Dynein heavy chain 14, axonemal (Axonemal beta dynein heavy chain 14) (Ciliary dynein heavy chain 14)
--------------------------	--

Product Name: Recombinant Human DNAH14 Protein
Catalog #: PEH9151



Note

For Research Use Only , Not for Diagnostic Use