


Product Name: Recombinant Human ANK3 Protein
Catalog #: PEH9153



Summary

Name	ANK3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q12955
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QMKPKLERII EVHIEKGNQA EPTEVIIRET KKHPEKEMYV YQKDLSRGDI NLKDFLPEKH DAFPCSEEQG QQEEEELTAE ESLPSYLESS RVNTPVSQEE DSRPSSAQLI SDDSYKTLKL LSQHSIEYHD DELSELRGES YRFAEKMLLS EKLDVSHSDT EESVTDHAGP PSELQGSDK RSREKIATAP KKEILSKIYK DVSENGVGKV SKDEHFDKVT VLHYSGNVSS PKHAMWMRFT EDRLDRGREK LIYEDRVDRT VKEAEKLTE VSQFFRDKTE KLNDELQSPE KKARPKNGKE YSSQSPTSSS

image



Background

Alternative Names Ankyrin-3 (ANK-3) (Ankyrin-G)

Product Name: Recombinant Human ANK3 Protein
Catalog #: PEH9153



Note

For Research Use Only , Not for Diagnostic Use