


Product Name: Recombinant Human CALCOCO2 Protein
Catalog #: PEH9156



Summary

Name	CACO2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13137
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRKDWIGIFR VGWKTTREYY TFMWVTLPID LNNKSAKQQE VQFKAYYLPK DDEYYQFCYV DEDGVVRGAS IPFQFRPENE EDILVVTQGEVEEIEQHMK ELCKENQELK DSCISLQKQN SDMQAELQKK QEELETLSI NKKLELVKE QKDYWETELL QLKEQNQKMS SENEKMGIRV DQLQAQLSTQ EKEMEKLQVQ DQDKTEQLEQ LKKENDHFL SLTEQRKDQK

image



Background

Alternative Names	Calcium-binding and coiled-coil domain-containing protein 2 (Antigen nuclear dot 52 kDa protein) (Nuclear domain 10 protein NDP52) (Nuclear domain 10 protein 52) (Nuclear dot protein 52)
--------------------------	--



Product Name: Recombinant Human CALCOCO2 Protein
Catalog #: PEH9156

Note

For Research Use Only , Not for Diagnostic Use