


Product Name: Recombinant Human GRAP Protein
Catalog #: PEH9159



Summary

Name	GRAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13588
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MESVALYSFQ ATESELA FN KGDTLKILNM EDDQNWYKAE LRGVEGFIPK NYIRVKPHPW YSGRISRQLA EEILMKRNHL GAFLIRESES SPGEFSVSVN YGDQVQHFKV LREASGKYFL WEEKFNSLNE LVDFYRTTTI AKKRQIFLRD EEPLLKSPGA CFAQAQFDFS AQDPSQLSFR RGDIIIEVLER PDPHWWRGRS CGRVGFFPRS YVQPVHL

image



Background

Alternative Names GRB2-related adapter protein

Note

Product Name: Recombinant Human GRAP Protein
Catalog #: PEH9159



For Research Use Only , Not for Diagnostic Use