


Product Name: Recombinant Human THOC5 Protein
Catalog #: PEH9164



Summary

Name	THOC5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13769
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKEITKCLEF KSKHEEIDLVSLEEFYKEAP PDISKAEVTMGDPHQQTLAR LDWELEQRKRLAEKYRECLSNKEKILKEIEVKKEYLSSLQ PRLNSIMQAS LPVQEYLFMPFDQAHKQYETARHLPPPLYVLFVQATAYGQACDKTLSVAI EGSVDEAKALFKPPEDSQDD ESDSDAEEEEQ

image



Background

Alternative Names	THO complex subunit 5 homolog (Functional spliceosome-associated protein 79) (fSAP79) (NF2/meningioma region protein pK1.3) (Placental protein 39.2) (PP39.2) (hTREX90)
--------------------------	---

Note

Product Name: Recombinant Human THOC5 Protein
Catalog #: PEH9164



For Research Use Only , Not for Diagnostic Use