


Product Name: Recombinant Human CAPS Protein
Catalog #: PEH9165



Summary

Name	CAYP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13938
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TTQLTQGPGAG RTLGQTQASC PEPRPSMDAV DATMEKLRAQ CLSRGASGIQ GLARFFRQLD RDGSRSLDAD EFRQGLAKLG LVLDQAEAEV VCRKWDRNGS GTLGLEEFRLR ALRPPMSQAR EAVIAAAFAK LDRSGDGVVT VDDLGRGVYSG RAHPKVRSGE WTEDEVLRRF LDNFDSSSEKD GQVTLAEFQD YYSGVVSASMN TDEEFVAMMT

image



Background

Alternative Names Calcyplosin (Calcyplosine)

Note

Product Name: Recombinant Human CAPS Protein
Catalog #: PEH9165



For Research Use Only , Not for Diagnostic Use