


Product Name: Recombinant Human TRIM29 Protein
Catalog #: PEH9167



Summary

Name	TRI29_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14134
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NGTKADGKDA KTTNGHGGEA AEGKSLGSAL KPGEGRSALF AGNEWRRPII QFVESGDDKN SNYFSMDSME GKRSPYAGLQ LGAACKPPVT FAEKGELRKS IFSESRKPTV SIMEPGETRR NSYPRADTGL FSRKSGSEE VLCDSCIGNK QKAVKSCLVC QASFCELHLK PHLEGAAFRD HQLLEPIRDF EARKCPVHGK TMELFCQTDQ TCICYLCMFQ EHKNHSTVTV EEAKAEKETE LSLQKEQLQL KIIEIEDEAE KWQKEKDRIK SFTTNEKAIL EQNFRDLVRD LEKQKEEVRA ALEQREQDAV DQVKVIMDAL DERAKVLHED

image



SDS

Background

Alternative Names	Tripartite motif-containing protein 29 (Ataxia telangiectasia group D-associated protein)
--------------------------	---

Product Name: Recombinant Human TRIM29 Protein
Catalog #: PEH9167



Note

For Research Use Only , Not for Diagnostic Use