


Product Name: Recombinant Human PPP2R5D Protein
Catalog #: PEH9170



Summary

Name	2A5D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14738
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PPKVAKCTAK PSSSGKDGGG ENTEEAQPQP QPQPQPQAQS QPPSSNKRPS NSTPPPTQLS KIKYSGGPQI VKKERRQSSS RFNLSKNREL QKLPALKDSP TQEREELFIQ KLRQCCVLFD FVSDPLSDLK FKEVKRAGLN EMVEYITHSR DVVTEAIYPE AVTMFSVNLF RTLPPSSNPT GAFFDPEEDE

image



Background

Alternative Names	Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit delta isoform (PP2A B subunit isoform B'-delta) (PP2A B subunit isoform B56-delta) (PP2A B subunit isoform PR61-delta) (PP2A B subunit isoform R5-delta)
--------------------------	--

Note

Product Name: Recombinant Human PPP2R5D Protein
Catalog #: PEH9170



For Research Use Only , Not for Diagnostic Use