


Product Name: Recombinant Human KIF19 Protein
Catalog #: PEH9180



Summary

Name	KIF19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2TAC6
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KIDEQTGRGQ ARGRQDRGDI RHIQAEVQLH SGQGEKAGMG QLREQLASAF QEQMDVRRRL LELENRAMEV QIDTSRHLLT IAGWKHEKSR RALKWREEQR KECYAKDDSE KSDTGDDQP DILEPPEVAA ARESIAALVD EQKQLRKQKL ALEQRCREL RRRLEETL PRRIGSEEQR EVLSLLCRVH ELEVENTEMQ SHALLRDGAL RHRHEAVRRL EQHRSLCDEI IQGQRQIIDD YNLAVPQRLE ELYEVYLREL EEGSLEQATI MDQVASRALQ

image



Background

Alternative Names Kinesin-like protein KIF19

Product Name: Recombinant Human KIF19 Protein
Catalog #: PEH9180



Note

For Research Use Only , Not for Diagnostic Use