


Product Name: Recombinant Human HKDC1 Protein
Catalog #: PEH9181



Summary

Name	HKDC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2TB90
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRV LALFQLT REQLVDVQAK MRAELEYGLK KKSHGLATVR MLPTYVCGLP DGTEKGKFLA LDLGGTNFRV LLVKIRSGRR SVRMYNKIFA IPLEIMQGTG EELFDHIVQC IADFLDYMGL KGASLPLGFT FSFPCRQMSI DKGTLLIGWTK GFKATDCEGE DVVDMLREAI KRRNEFDLDI VAVVNDTVGT MMTCGYEDPN CEIGLIAGTG SNMCYMEDMR NIEMVEGGEG KMCINTEWGG FGDNGCIDDI WTRYDTEVDE GSLNPGKQRY EKMTSGMYLG EIVRQILIDL TKQGLLFRGQ ISERLRTRGI FETKFLSQIE SDR LALLQVR RILQQGLGLDS TCEDSIVVKE VCGAVSRRAA QLCGAGLAAI VEKRREDQGL EHLRITVGVD GTLYKLHPHF SRILQETVKE LAPRCDVTFM LSEDGSGKGA ALITAVAKRL

image



Background

Product Name: Recombinant Human HKDC1 Protein
Catalog #: PEH9181



Alternative Names

Hexokinase HKDC1 (EC 2.7.1.1) (Hexokinase domain-containing protein 1)

Note

For Research Use Only , Not for Diagnostic Use