


Product Name: Recombinant Human CTU2 Protein
Catalog #: PEH9182



Summary

Name	CTU2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2VPK5
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RAGDAFCRDC FKA FYVHKFR AMLGKNRLIF PGEKVLLAWS GGPSSSSMVW QVLEGLSQDS AKRLRFVAGV IFVDEGAACG QSLEERSKTL AEVKPILQAT GFPWHVVALE EVFSLPPSVL WCSAQELVGS EGAYKAAVDS FLQQQHVLGA GGGPGPTQGE EQPPQPPLDP QNLARPPAPA QTEALSQ LFC SVRTLTAKEE LLQTLRTHLI LHMARAHGYS KVMTGDSCTR LAIKLMTNLA LGRGAFLAWD TGFSDERHGD

image



Background

Alternative Names Cytoplasmic tRNA 2-thiolation protein 2 (Cytosolic thiouridylase subunit 2)

Product Name: Recombinant Human CTU2 Protein
Catalog #: PEH9182



Note

For Research Use Only , Not for Diagnostic Use