


Product Name: Recombinant Human PUS10 Protein
Catalog #: PEH9183



Summary

Name	PUS10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3MIT2
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLETEKDELI LEVMNPPPKK IRLQELEDSI DNLSQNGEGR ISVSHVGSTA SKNSNLNVCN VCLGILQEFC

image



Background

Alternative Names	tRNA pseudouridine synthase Pus10 (Hup10) (EC 5.4.99.25) (Coiled-coil domain-containing protein 139) (tRNA pseudouridine 55 synthase) (Psi55 synthase) (tRNA pseudouridylate synthase) (tRNA-uridine isomerase)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use