


Product Name: Recombinant Human RABL6 Protein
Catalog #: PEH9184



Summary

Name	RABL6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q3YEC7
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDQAPGRDKN IPAGLQSMNQ ALQRRFAKGV QYNMKIVIRG DRNTGKTALW HRLQGRPFVE EYIPTQEIQV TSIHWSYKTT DDIVKVEVWD VVDKGKCKKR GDGLKMENDP QEAESMALD AEFLDVYKNC NGVMMFDIT KQWTFNYILR ELPKVPTHVP VCVLGNYRDM GEHRVILPDD VRDFIDNLD RPPGSSYFRYA ESSMKNSFGL KYLHKFFNIP FLQLQRETL RQLETNQLDM DATLEELSVQ QETEDQNYGI FLEMMEARSR GHASPLAANG QSPSPGSQSP VVPAGAVSTG SSSPGTPQPA PQLPLNAAPP

image



Background

Alternative Names	Rab-like protein 6 (GTP-binding protein Parf) (Partner of ARF) (Rab-like protein 1) (RBEL1)
--------------------------	---

Product Name: Recombinant Human RABL6 Protein
Catalog #: PEH9184



Note

For Research Use Only , Not for Diagnostic Use