


Product Name: Recombinant Human FILIP1L Protein
Catalog #: PEH9185



Summary

Name	FIL1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q4L180
Host	E.coli
Species	Human
Molecular Mass	63 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EFICLLEQEC ERLKKLIDQE IKSQEEKEQE KEKRVTTLKE ELTKLKSFAL MNVDEQQRLT AQLTLQRQKI QELTTNAKET HTKLALAEAR VQEEEQKATR LEKELQTQTT KFHQDQDTIM AKLTNEDSQN RQLQQKLAAL SRQIDELEET NRSLRKAEED LQDIKEKISK GEYGNAGIMA EVEELRKRVL DMEGKDEELI KMEEQCRDLN KRLERETLQS KDFKLEVEKL SKRIMALEKL EDAFNKSKQE CYSLKCNLEK ERMTTKQLSQ ELESKVRIK ELEAIESRLE KTEFTLKEDL TKLKTTLVMF VDERKTMSEK LKKTEDKLQA ASSQLQVEQN KVTTVTEKLI EETKRALKSK TDVEEKMYSV TKERDDLKNK LKAEEEKGND LLSRVNMLKN RLQSLEAIEK DFLKNKLNQD SGKSTTALHQ ENNIKIKELSQ EVERLKLKLK DMKAIEDDLM KTEDEYETLE RRYANERDKA QFLSKELEHV KMELAKYKLA EKTETSHEQW LFKRLQEEEA

image



Background

Product Name: Recombinant Human FILIP1L Protein
Catalog #: PEH9185



Alternative Names

Filamin A-interacting protein 1-like (130 kDa GPBP-interacting protein) (90 kDa GPBP-interacting protein) (Protein down-regulated in ovarian cancer 1) (DOC-1)

Note

For Research Use Only , Not for Diagnostic Use