


Product Name: Recombinant Human HSP90AB4P Protein
Catalog #: PEH9190

Summary

Name	H90B4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q58FF6
Host	E.coli
Species	Human
Molecular Mass	59 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKEIFLQELI SNASDALDKI RYESLTDP SK LDGGKELKID IIPNPRECIL TLVNTGIGMT KADLINNLGA IAKSGTEAFM EAFQSCAEIS MIGQFGVGFY SAYLVAEKVA ITKHNDDEEQY SWVSSAGSSF TLHVDHGEP I DRDTKVILHL KEDQTEYLEE RWVKEVVKKH PQFIGCLIAV YLEKEPEKEI SDDEEEKGEK EEEDKDDKEK PKTEDVGSDE EDDTDKNNKK KTKKIKEYT DREELNQTKP IWTRNPDDIT QEECGEFYKS LTSAWEDHLA VKQFPVEEQE ENEQLCVHHV WIMDSFDDLM PEYVGFVRED KENNKKLDEV FSKISWLG I H EDSINWRHLS ELLWSHTFQS GDEMTSLSEY VSCMKEAQKS ICDIIGECKE QVANS AFVEQ EWKKGFEVIY MSEPIDEYCV QQLKEFDGKS LLSVTKEGLE LPEDEEEKKI MEESNVKFEN LCRLMKEILD KKVERTISS RLVSSPCRIV

image



Background



Product Name: Recombinant Human HSP90AB4P Protein
Catalog #: PEH9190

Alternative Names

Putative heat shock protein HSP 90-beta 4

Note

For Research Use Only , Not for Diagnostic Use