


Product Name: Recombinant Human HSP90AA5P Protein
Catalog #: PEH9191



Summary

Name	HS905_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q58FG0
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KRKKKKESND KPEIEDVGSD EEEKKDADK KKKKSKEYI DQELNKTPI WTRNPDAITN EEEYGEFHQSL TNNWEDHLAV KHFSVEGQLE ELKDSRRVMK ANQKHIIYIT GETKDQVANS AFVECLQKHG LEVIYMIELI DKYCVQQLKE LESKTVVSVA KEGLELPEDE EEKKKQEEKK TKFENLCKIM

image



Background

Alternative Names	Putative heat shock protein HSP 90-alpha A5 (Heat shock protein 90-alpha E) (Heat shock protein 90Ae)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HSP90AA5P Protein
Catalog #: PEH9191

