


Product Name: Recombinant Human BEX5 Protein
Catalog #: PEH9193



Summary

Name	BEX5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5H9J7
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VEKAPVQNEA PALGGGEYQE PGGNVKGVWA PPAPGFGE DV PNRLVDNIDM IDGDGDDMER FMEEMRELRR KIRELQLRYS LRILIGDPPH HDHHDEFCLM

image



Background

Alternative Names Protein BEX5 (Brain-expressed X-linked protein 5) (NGFRAP1-like protein 1) (Nerve growth factor receptor-associated protein 2)

Note

For Research Use Only , Not for Diagnostic Use