


Product Name: Recombinant Human TAF7L Protein
Catalog #: PEH9194



Summary

Name	TAF7L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5H9L4
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GREEKCVWKH GITPLKNVR KKRFRKTQKK VPDVKEMEKS SFTEYIESPD VENEVKRLLR SDAEAVSTRW EVIAEDGTKE IESQGSIPGF LISSGMSSHK QGHTSSEYDM LREMFSDSRS NNDDDEDEDD EDEDEDEDED EDEDKEEEEE DCSEEYLERQ LQAEFIESGQ YRANEGTSSI VMEIQKQIEK KEKKLHKIQN KAQRQKDLIM KVENLTLKNH FQSVLEQLEL

image



Background

Alternative Names	Transcription initiation factor TFIID subunit 7-like (Cancer/testis antigen 40) (CT40) (RNA polymerase II TBP-associated factor subunit Q) (TATA box-binding protein-associated factor 50 kDa) (Transcription initiation factor TFIID 50 kDa subunit)
--------------------------	---

Product Name: Recombinant Human TAF7L Protein
Catalog #: PEH9194



Note

For Research Use Only , Not for Diagnostic Use