


Product Name: Recombinant Human EEF1AKMT2 Protein
Catalog #: PEH9197



Summary

Name	EFMT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5JPI9
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAVAARSDKG SPGEDGFVPS ALGTREHWDA VYERELQTFR EYGDTGEIWF GEESMNRLIR WMQKHKIPLD ASVLDIGTGN GVFLVELAKF GFSNITGIDY SPSAIQLSGS IIEKEGLSNI KLVEDFLNL STQLSGFHIC IDKGTFDAIS LNPDNAIEKR KQYVKSLSRV LKVKGFFLIT SCNWTKEELL NEFSEGWSTV AGFWLTAALT

image



Background

Alternative Names	EEF1A lysine methyltransferase 2 (EC 2.1.1.-) (Methyltransferase-like protein 10) (Protein-lysine N-methyltransferase METTL10)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human EEF1AKMT2 Protein
Catalog #: PEH9197

