


Product Name: Recombinant Human SAMD9 Protein
Catalog #: PEH9199



Summary

Name	SAMD9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5K651
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSKNAPKDQT VSQKERRETS KQKQKGKENP DMANPSAMST TAKGSKSLKV ELIEDKIDYT KERQPSIDLT CVSYPFDEFS NPYRYKLDFS LQPETGPGNL IDPIHEFKAF TNTATATEED VKMKFSNEVF RFASACMNSR TNGTIHFGVK DKPHGKIVGI

image



Background

Alternative Names	Sterile alpha motif domain-containing protein 9 (SAM domain-containing protein 9)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SAMD9 Protein
Catalog #: PEH9199

