


Product Name: Recombinant Human PERM1 Protein
Catalog #: PEH9201



Summary

Name	PERM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5SV97
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPRRAKLHTV STTVWEALPD VSRAKSDMAV STPASEPQPD RDMAVSTPAS EPQSDRDMAV STPASEPQPD

image



Background

Alternative Names	PGC-1 and ERR-induced regulator in muscle protein 1 (PPARGC1 and ESRR-induced regulator in muscle 1) (Peroxisome proliferator-activated receptor gamma coactivator 1 and estrogen-related receptor-induced regulator in muscle 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PERM1 Protein
Catalog #: PEH9201

