


Product Name: Recombinant Human DCAF12 Protein
Catalog #: PEH9204



Summary

Name	DCA12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T6F0
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QITKIPLKD REPGGVTTQQG CGIHAIELNP SRTLLATGGD NPNSLAIRL PTLDPVCGVD DGHKDWIFS AWISDTMAVS GSRDGSMGLW EVTDDVLTKS DARHNVS RVP VYAHITHKAL KDIPKEDTNP DNCKVRALAF NNKNKELGAV

image



Background

Alternative Names	DDB1- and CUL4-associated factor 12 (Centrosome-related protein TCC52) (Testis cancer centrosome-related protein) (WD repeat-containing protein 40A)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human DCAF12 Protein
Catalog #: PEH9204

