



Product Name: Recombinant Human RNF187 Protein
Catalog #: PEH9206



Summary

Name	RN187_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TA31
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDLNDARDLH GQAESAAVW KGHVMDRRKK ALTDYKKLRA FFVEEEHFL QEAEKEEGLP EDELADPTER FRSLQAVSE

image



Background

Alternative Names	E3 ubiquitin-protein ligase RNF187 (EC 2.3.2.27) (RING domain AP1 coactivator 1) (RACO-1) (RING finger protein 187) (RING-type E3 ubiquitin transferase RNF187)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use