


Product Name: Recombinant Human NBPF12 Protein
Catalog #: PEH9207



Summary

Name	NBPFC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TAG4
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	INEKLRPQLA EKKQQFRSLK EKCFVTQLAG FLAKQQNKYK YEECKDLIKS MLRNEQLQFKE EKLAELQKQA EELRQYKVLV HSQERELTQL REKLREGRDA SRSLNEHLQA LLTPDEPDKS QGQDLQEQLA EGCRLAQHLV QKLSPENDE EDEDVQVEED EKVLESSPR EMQKAEESKV PEDSLEECAI TCSNSHGPCD SNQPHKNIKI TFEEDKVNSS LVVDRESSHD ECQDALNILP VPGPTSSATN VSMVVSAGPL SSEKAEMNIL EINEKLRPQL AEKKQQFRSL KEKCFVTQVA CFLAKQQNKY KYEICKDLLK SMLRNEQLQFK EEKLAELQKQ AEELRQYKVL VHSQERELTQ LREKLREGRD ASRSLNEHLQ ALLTPDEPDK SQGQDLQEQL AEGCRLAQHL

image



Background

Product Name: Recombinant Human NBPF12 Protein
Catalog #: PEH9207



Alternative Names

Neuroblastoma breakpoint family member 12 (Chromosome 1 amplified sequence 1)

Note

For Research Use Only , Not for Diagnostic Use