


Product Name: Recombinant Human DCAF8 Protein
Catalog #: PEH9208



Summary

Name	DCAF8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TAQ9
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGRETSSGIE VEASDLISL TGDDGGPNRT STESRGTDTE SSGEDKSDS MEDTGHYSIN DENRVHDRSE EEEEEEEEEEE EQPRRRVQR KRANRDQDSS DDERALEDWV SSETSALPRP RWQALPALRE RELGSSARFV YEACGARVFV QRFRLQHGLE GHTGCVNTLH FNQRGTWLAS GSDDLKVVVW DWVRRQPVLD FESGHKSNVF QAKFLPNSGD STLAMCARDG QVRVAELSAT QCCKNTKRVA QHKGASHKLA LEPDSPCTFL SAGEDAVVFT IDLRQDRPAS KLVVTKEKEK

image



Background

Alternative Names DDB1- and CUL4-associated factor 8 (WD repeat-containing protein 42A)

Product Name: Recombinant Human DCAF8 Protein
Catalog #: PEH9208



Note

For Research Use Only , Not for Diagnostic Use