


Product Name: Recombinant Human STPG1 Protein
Catalog #: PEH9210



Summary

Name	STPG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TH74
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TGKHPRRASE VQKGFTAAYP TQSSIPFKSQ ASVIPSEKK GFNSQAKRFP HKKNDIPGPG FYNVIHQSPV

image



Background

Alternative Names O(6)-methylguanine-induced apoptosis 2 (MAPO2) (Sperm-tail PG-rich repeat-containing protein 1)

Note

For Research Use Only , Not for Diagnostic Use