



Product Name: Recombinant Human MEX3C Protein
Catalog #: PEH9212



Summary

Name	MEX3C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5U5Q3
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DEEEGEEAEL DGDLEEEEL EEAEEDRSS LLLLSPPAAT ASQTQQIPGG SLGSVLLPAA RFDAREAAAA AAAAGVLYGG DDAQGMMMAAM LSHAYGPGGC GAAAAALNGE QAALLRRKSV NTTECVVPS SEHVAEIVGR QGCKIKALRA KTNTYIKTPV RGEPIFVVT GRKEDVAMAK REILSAAEHF SMIRASRNKN

image



Background

Alternative Names	RNA-binding E3 ubiquitin-protein ligase MEX3C (EC 2.3.2.27) (RING finger and KH domain-containing protein 2) (RING finger protein 194) (RING-type E3 ubiquitin transferase MEX3C)
--------------------------	---

Note

Product Name: Recombinant Human MEX3C Protein
Catalog #: PEH9212



For Research Use Only , Not for Diagnostic Use