


Product Name: Recombinant Human MLIP Protein
Catalog #: PEH9215



Summary

Name	MLIP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VWP3
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EESSDKSPET VNRSKSNLYL TLNAGSQQER DQAKLTCPSE VSGTILQERE FEANKLQGMQ QSDLFKAIEYV LIVDSEGEDE AASRKVEQGP PGGIGTAAVR PKSLAIISSL VSDVVRPKTQ GTDLKTSSHP EMLHGMAPQQ KHGQQYKTKS SYKAFAAIPT NTLLEQKAL DEPAKTESVS KDNTLEPPVE LYFPAQLRQQ TEELCATIDK VLQDSLMSHS SDSPSRSPKT LLGSDTVKTP TTLPRAAGRE

image



Background

Alternative Names	Muscular LMNA-interacting protein (Cardiac Isl1-interacting protein) (CIP) (Muscular-enriched A-type laminin-interacting protein)
--------------------------	--

Note

Product Name: Recombinant Human MLIP Protein
Catalog #: PEH9215



For Research Use Only , Not for Diagnostic Use