


Product Name: Recombinant Human STYXL2 Protein
Catalog #: PEH9216



Summary

Name	STYL2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VZP5
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEPIHLSSAI AAKQIINEEL KPPGVRADAE CPGMLESAEQ LLVEDLYNRV REKMDDTSLY NTPCVLDLQR ALVQDRQEAP WNEVDEVWPN VFIAEKSVAV NKGRLKRLGI THILNAAHGT GVYTGPEFYT GLEIQYLGVE VDDFPEVDIS QHFRKASEFL DEALLTYRGK VLVSSEMGIS RSAVLVVAYL MIFHNMAILE ALMTVRKKRA IYPNEGFLKQ LRELNEKLME EREEDYGREG GSAEAEEGEG TGSM LGARVH ALTVEEEDDS ASHLSGSSLG KATQASKPLT LIDEEEEKL YEQWKKGQGL LSDKVPQDGG GWRSSASSGQG GEELEDEDVE RIIQEWQSRN

image



Background

Alternative Names	Serine/threonine/tyrosine-interacting-like protein 2 (Inactive dual specificity phosphatase 27)
--------------------------	---

Product Name: Recombinant Human STYXL2 Protein
Catalog #: PEH9216



Note

For Research Use Only , Not for Diagnostic Use