


Product Name: Recombinant Human BLTP3A Protein
Catalog #: PEH9219



Summary

Name	URFB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6BDS2
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DHYQYLALLR LKEVLQRLQE QLTKDTESMT GSPLQNQTAC IGVLFPSAEV ALLMHPAPGA VDADSAGSDS TSLVDSELSF SEDRELKSDA SSDQGPASPE KVLLEESSIEN QDVSQERPHS

image



Background

Alternative Names	UHRF1-binding protein 1 (ICBP90-binding protein 1) (Ubiquitin-like containing PHD and RING finger domains 1-binding protein 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use