



## Summary

<b>Name</b>	PPR18_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6NYC8
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	AIGPVHQNRF IRQERQQQQ QQQRSEELLA ERKPGPLEAR ERRPSPGEMR DQSPKGRESR EERLSPRETR ERLGIGGAQ ELSLRPLEAR DWRQSPGEVG DRSSLSEAW KWRLSPGETP ERSRLAESR EQSPRRKEVE SRLSPGESAY QKLGLTEAHK WRPDSRESQE QSLVQLEATE WRLRSGEERQ DYSEECGRKE EWPVPGVAPK ETAELSETLT REAQGNSSAG VEAAEQRPVE DGERGMKPT GWKWTLSNGK AREWTPRDIE AQTQKPEPPE SAEKLLLESPG VEAGEGAEAK EEAGAQRPL

## image



## Background

<b>Alternative Names</b>	Phostensin (Protein phosphatase 1 F-actin cytoskeleton-targeting subunit) (Protein phosphatase 1 regulatory subunit 18)
--------------------------	--

**Product Name: Recombinant Human PPP1R18 Protein**  
**Catalog #: PEH9221**



---

**Note**

For Research Use Only , Not for Diagnostic Use