


**Product Name: Recombinant Human CCDC62 Protein**  
**Catalog #: PEH9223**



## Summary

|                              |   |
|------------------------------|---|
| <b>Name</b>                  | CCD62_HUMAN   |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE   |
| <b>Endotoxin level</b>       | Please contact with the lab for this information  |
| <b>Tag</b>                   | C-His   |
| <b>Accession #</b>           | Q6P9F0  |
| <b>Host</b>                  | E.coli  |
| <b>Species</b>               | Human   |
| <b>Molecular Mass</b>        | 31 kDa  |
| <b>Buffer</b>                | Phosphate buffered saline   |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.  |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.  |
| <b>Sequence</b>              | LEARNETLSN TLVELSAQVG QLQAREQALT TMIKDKDI IEAVNHIADC<br>SGKFKMLEHA LRDAKMAETC IVKEKQDYKQ KLKALKIEVN KLKEDLNEKT<br>TENNEQREEI IRLKQEKSL HDPELLFTVER EKRKDELLNI AKSKQERTNS<br>ELHNLRQIYV KQQSDLQFLN FNVENSQELI QMYDSKMEES KALDSSRDMC<br>LSDLENNHPK VDIKREKNQK SLFKDQKFEA MLVQQNRSK SSCDECKEKK |

## image



## Background

**Alternative Names** Coiled-coil domain-containing protein 62 (Protein TSP-NY) (Protein aaa)

## Note

**Product Name: Recombinant Human CCDC62 Protein**  
**Catalog #: PEH9223**



---

For Research Use Only , Not for Diagnostic Use