


Product Name: Recombinant Human BLOC1S3 Protein
Catalog #: PEH9225



Summary

Name	BL1S3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6QNY0
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRPETVVPG EATETDSERS ASSSEEEELY LGPSGPTRGR PTGLRVAGEA AETDSEPEPE PEPTAAPRDL PPLVVQRESA EEAWGTEEAP APAPARSLQ LRLAESQARL DHDVAAAVSG

image



Background

Alternative Names	Biogenesis of lysosome-related organelles complex 1 subunit 3 (BLOC-1 subunit 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use