


Product Name: Recombinant Human STAC2 Protein
Catalog #: PEH9228



Summary

Name	STAC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZMT1
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDGEQSIRSS EEGPGDSASP VFTAPAESEG PGPEEKSPGQ QLPKATLRKD VGPMYSYVAL YKFLPQENND LALQPGDRIM LVDDSNEDWW KGKIGDRVGF FPANFVQRVR PGENVWRCCQ PFSGNKEQGY MSLKENQICV GVGRSKDADG FIRVSSGKKR GLVPVDALTE

image



Background

Alternative Names SH3 and cysteine-rich domain-containing protein 2 (24b2/STAC2) (Src homology 3 and cysteine-rich domain-containing protein 2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human STAC2 Protein
Catalog #: PEH9228

