


Product Name: Recombinant Human CEP128 Protein
Catalog #: PEH9230



Summary

Name	CE128_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZU80
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERALEKQSET VDELTGKNNQ ILKEKDELKT QLYAALQQIE NLRKELNDVL TKRALQEEEL HSKEEKL RDI KSHQADLELE VKNSLDTIHR LESELKKQSK IQSQMKVEKA HLEEEIAELK KSQAQDKAKL LEMQESIKDL SAIRADLANK LAEERAKKA

image



Background

Alternative Names Centrosomal protein of 128 kDa (Cep128)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CEP128 Protein
Catalog #: PEH9230

