


**Product Name: Recombinant Human ARHGAP27 Protein**  
**Catalog #: PEH9231**



## Summary

<b>Name</b>	RHG27_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q6ZUM4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	35 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	EMQPGLSPGS PGDPRPPTPE TDYPESLTSY PEEDYSPVGS FGEPGPTSPL TTPPGWSCHV SQDKQMLYTN HFTQEQQWVRL EDPHGKPYFY NPEDSSVRWE LPQVPVPAPR SIHKSSQDGD TPAQASPPEE KVPAELDEVG SWEEVSPATA AVRTKTLDKA GVLHRTKTAD KGKRLRKKHW SASWTVLEGG VLTFFKDSKT SAAGGLRQPS KFSTPEYTVE LRGATLSWAP KDKSSRKNVL ELRSRDGSEY LIQHDSEAI STWHKAI AQG IQELSAELPP EESESRVDF GSSERLGSWQ

## image



## Background

<b>Alternative Names</b>	Rho GTPase-activating protein 27 (CIN85-associated multi-domain-containing Rho GTPase-activating protein 1) (Rho-type GTPase-activating protein 27) (SH3 domain-containing protein 20)
--------------------------	--

**Product Name: Recombinant Human ARHGAP27 Protein**  
**Catalog #: PEH9231**

---



**Note**

For Research Use Only , Not for Diagnostic Use