



Product Name: Recombinant Human MICAL2 Protein
Catalog #: PEH9232



Summary

Name	MICLK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZW33
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RDLFGSPKRK VLPEDSAQAL EKLLQPFKST SLRQAAPPPP PTPPPPPPPP TAGGADSKNF PLRAQVTEAS SSASSTSSSS ADEEFDPQLS LQLKEKKTLR RRKKLEKAMK QLVKQEELKR LYKAQAIQRQ LEEVEERQRA SEIQGVRLEK ALRGEADSGT QDEAQLLQEW FKLVLKKNKL MRYSELLIM AQELELEDHQ SRLEQKLREK MLKEESQKDE KDLNEEQEVF TELMQVIEQR DKLVDLSLEEQ RIREKAEDQH

image



Background

Alternative Names	MICAL C-terminal-like protein (Mical-cL) (ERK2-binding testicular protein 1) (Ebitein-1)
--------------------------	--

Product Name: Recombinant Human MICAL2 Protein
Catalog #: PEH9232



Note

For Research Use Only , Not for Diagnostic Use