



Summary

Name	ANKS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZW76
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QYGHMKIVAL MDTYSPSLPK SLYRSPEKYE DLSSDESCP APQRQRPCRK KGVSIHEGPR ALARITGIGL GGRAPRPRYE QAPPRGYVTF NSSGENPLEE EGLCCRDVTS PINERDVES SSSSREEHA FCANLGPVQS SSSSEGLARA QGLSSEASVE SNEDSDHACK SSARKQAKSY MKTKNPDSQW PPRAATDREG FLAESSPQTQ

image



Background

Alternative Names Ankyrin repeat and SAM domain-containing protein 3

Note

Product Name: Recombinant Human ANKS3 Protein
Catalog #: PEH9233



For Research Use Only , Not for Diagnostic Use